

Specification Sheet



Product Name: Lesinurad

Item No.: MED10454

Purity Specification: >99%

Molecular Formula : C₁₇H₁₃BRN₃NAO₂S

CASNumber: 878672-00-5

FormulaWeight: 426.263


Storage: -20

Test

Specification

	Min.	Max.
PURITY	99%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816