

Specification Sheet



Product Name: Pexidartinib (PLX3397)

Item No.: MED10688

Purity Specification: >98%


Molecular Formula : C₂₀H₁₅CLF₃N₅

CASNumber: 1029044-16-3 **FormulaWeight:** 417.815

Storage: -20

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816