

Specification Sheet



Product Name: Bax inhibitor peptide P5

Item No.: MED11236

Purity Specification: 98%


Molecular Formula : C₂₇H₄₈N₆O₈S

CASNumber: 579492-83-4 **FormulaWeight:** 616.770

Storage: Desiccate at -20 °C

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816