

# Specification Sheet



**Product Name:** Dynamin inhibitory peptide

**Item No.:** MED11370

**Purity Specification:** 98%

**Molecular Formula :** C<sub>47</sub>H<sub>80</sub>N<sub>18</sub>O<sub>14</sub>

**CAS Number:** 251634-21-6      **Formula Weight:** 1121.250


**Storage:** Desiccate at -20 °C

## Test

## Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: [www.med-bio.cn](http://www.med-bio.cn) Tel: 400-186-1816