

# Specification Sheet



**Product Name:** Iniparib (BSI-201)

*Item No.:* MED11676

**Purity Specification:** 98%


**Molecular Formula :** C<sub>7</sub>H<sub>5</sub>IN<sub>2</sub>O<sub>3</sub>

**CASNumber:** 160003-66-7    **FormulaWeight:** 292.031

**Storage:** -20 ° C

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn    Tel:400-186-1816