

Specification Sheet



Product Name: LMK 235

Item No.: MED11971

Purity Specification: 98%


Molecular Formula : C₁₅H₂₂N₂O₄

CAS Number: 1418033-25-6 **Formula Weight:** 294.346

Storage: -20 ° C

| Test | Specification | |
|--------|---------------|------|
| | Min. | Max. |
| PURITY | 98% | 100% |

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816