

Specification Sheet



Product Name: Clonidine HCl

Item No.: MED12581

Purity Specification: 98%

Molecular Formula : C₉H₁₀CL₃N₃


CASNumber: 4205-91-8

FormulaWeight: 266.555

Storage: -20 ° C

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816