

Specification Sheet



Product Name: Pancreatic Polypeptide (human)

Item No.: MED12698

Purity Specification: 98%

Molecular Formula : C₁₈₅H₂₈₇N₅₃O₅₄S₂

CAS Number: 75976-10-2

Formula Weight: 4181.711


Storage: Desiccate at -20 °C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816