

Specification Sheet



Product Name: Apremilast (CC-10004)

Item No.: MED14759

Purity Specification: 98%

Molecular Formula : C₂₂H₂₄N₂O₇S


CASNumber: 608141-41-9

FormulaWeight: 460.500

Storage: -20 ° C

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816