

Specification Sheet



Product Name: Pneumocandin B0

Item No.: MED15304

Purity Specification: 98%

Molecular Formula : C₄₉H₇₈N₈O₁₇

CASNumber: 135575-42-7

FormulaWeight: 1051.187


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816