

Specification Sheet



Product Name: Isochlorogenic acid C

Item No.: MED15919

Purity Specification: 98%

Molecular Formula : C₂₇H₂₈O₁₀

CAS Number: 32451-88-0

Formula Weight: 512.505


Storage: -20 ° C

Test

Specification

| | Min. | Max. |
|--------|------|------|
| PURITY | 98% | 100% |

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816