

Specification Sheet



Product Name: shikonofuran A

Item No.: MED16040

Purity Specification: 98%

Molecular Formula : C₁₈H₂₀O₅


CASNumber: 85022-66-8

FormulaWeight: 316.348

Storage: -20 ° C

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816