

Specification Sheet



Product Name: Prosapogenin CP6

Item No.: MED16101

Purity Specification: 98%

Molecular Formula : C₄₆H₇₄O₁₆

CASNumber: 72629-76-6

FormulaWeight: 883.070


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816