

Specification Sheet



Product Name: Sanggenone C

Item No.: MED16388

Purity Specification: 98%

Molecular Formula : C₄₀H₃₆O₁₂

CAS Number: 80651-76-9

Formula Weight: 708.707


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816