

Specification Sheet



Product Name: Hypericin

Item No.: MED16430

Purity Specification: 98%

Molecular Formula : C₃₀H₁₆O₈

CAS Number: 548-04-9

Formula Weight: 504.443


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816