

Specification Sheet



Product Name: Sophoridine

Item No.: MED16448

Purity Specification: 98%

Molecular Formula : C₁₅H₂₄N₂O


CASNumber: 6882-68-4

FormulaWeight: 248.364

Storage: -20 ° C

| Test | Specification | |
|--------|---------------|------|
| | Min. | Max. |
| PURITY | 98% | 100% |

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816