

Specification Sheet



Product Name: PTIQ

Item No.: MED17848

Purity Specification: 98%


Molecular Formula : C₁₃H₁₇NO₃

CASNumber: 1032822-42-6 **FormulaWeight:** 235.279

Storage: Store at RT

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816