

# Specification Sheet



**Product Name:** L-Biopterin

*Item No.:* MED17963

**Purity Specification:** 98%

**Molecular Formula :** C<sub>9</sub>H<sub>11</sub>N<sub>5</sub>O<sub>3</sub>

**CASNumber:** 22150-76-1

**FormulaWeight:** 237.215


**Storage:** -20 ° C

## Test

## Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816