

# Specification Sheet



**Product Name:** TG100713

*Item No.:* MED18161

**Purity Specification:** 98%

**Molecular Formula :** C<sub>12</sub>H<sub>10</sub>N<sub>6</sub>O


**CASNumber:** 925705-73-3

**FormulaWeight:** 254.247

**Storage:** -20 ° C

| Test   | Specification |      |
|--------|---------------|------|
|        | Min.          | Max. |
| PURITY | 98%           | 100% |

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816