

Specification Sheet



Product Name: TNP-470

Item No.: MED18550

Purity Specification: 98%


Molecular Formula : C₁₉H₂₈CLNO₆

CASNumber: 129298-91-5 **FormulaWeight:** 401.882

Storage: -20 ° C

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816