

# Specification Sheet



**Product Name:** TG003

*Item No.:* MED19784

**Purity Specification:** 98%

**Molecular Formula :** C<sub>13</sub>H<sub>15</sub>NO<sub>2</sub>S

**CAS Number:** 300801-52-9

**Formula Weight:** 249.329


**Storage:** -20 ° C

## Test

## Specification

	<b>Min.</b>	<b>Max.</b>
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: [www.med-bio.cn](http://www.med-bio.cn) Tel: 400-186-1816