

Specification Sheet



Product Name: Heparin

Item No.: MED20148

Purity Specification: 98%

Molecular Formula : C₂₆H₄₂N₂O₃₇S₅

CASNumber: 9005-49-6

FormulaWeight: 1134.928


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816