

Specification Sheet



Product Name: Icaritin

Item No.: MED21516

Purity Specification: 98%


Molecular Formula : C₂₁H₂₀O₆

CAS Number: 118525-40-9 **Formula Weight:** 368.380

Storage: -20 ° C

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816