

# Specification Sheet



**Product Name:** DL-Serine

*Item No.:* MED22317

**Purity Specification:** 98%

**Molecular Formula :** C<sub>3</sub>H<sub>7</sub>NO<sub>3</sub>


**CASNumber:** 302-84-1

**FormulaWeight:** 105.093

**Storage:** -20 ° C

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: [www.med-bio.cn](http://www.med-bio.cn) Tel: 400-186-1816