

Specification Sheet



Product Name: Isosorbide 5-mononitrate

Item No.: MED22487

Purity Specification: 98%

Molecular Formula : C₆H₇NO₆


CASNumber: 16051-77-7

FormulaWeight: 189.123

Storage: -20 ° C

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816