

Specification Sheet



Product Name: 3-Indolepropionic acid

Item No.: MED22684

Purity Specification: 98%

Molecular Formula : C₁₁H₁₁NO₂

CAS Number: 830-96-6

Formula Weight: 189.211


Storage: -20 ° C

Test

Specification

	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816