

Specification Sheet



Product Name: Prostaglandin E3

Item No.: MED23884

Purity Specification: 98%

Molecular Formula : C₂₀H₃₀O₅

CASNumber: 802-31-3

FormulaWeight: 350.449


Storage: -20

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816