

Specification Sheet



Product Name: LCK Inhibitor

Item No.: MED23901

Purity Specification: 98%

Molecular Formula : C₂₃H₂₂N₄O

CASNumber: 213743-31-8

FormulaWeight: 370.447


Storage: -20

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816