

Specification Sheet



Product Name: Losartan-d4

Item No.: MED24314

Purity Specification: 98%


Molecular Formula : C₂₂H₂₃CLN₆O

CASNumber: 1030937-27-9 **FormulaWeight:** 426.935

Storage: 2-8

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816