

Specification Sheet



Product Name: D-Kynurenine

Item No.: MED26617

Purity Specification: 98%

Molecular Formula : C₁₀H₁₂N₂O₃


CASNumber: 13441-51-5

FormulaWeight: 208.214

Storage: -20

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816