

Specification Sheet



Product Name: Taleranol

Item No.: MED26633

Purity Specification: 98%

Molecular Formula : C₁₈H₂₆O₅

CASNumber: 42422-68-4

FormulaWeight: 322.396


Storage: -20

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816