

Specification Sheet



Product Name: A2AR Agonist-1

Item No.: MED26659

Purity Specification: 98%

Molecular Formula : C₂₀H₂₂N₆O₄

CASNumber: 41552-95-8

FormulaWeight: 410.426


Storage: -20

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816