

Specification Sheet



Product Name: Paquinimod

Item No.: MED26674

Purity Specification: 98%

Molecular Formula : C₂₁H₂₂N₂O₃

CASNumber: 248282-01-1

FormulaWeight: 350.411


Storage: -20

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816