

Specification Sheet



Product Name: Pleconaril

Item No.: MED26687

Purity Specification: 98%

Molecular Formula : C₁₈H₁₈F₃N₃O₃

CASNumber: 153168-05-9

FormulaWeight: 381.349


Storage: -20

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816