

Specification Sheet



Product Name: Hexarelin

Item No.: MED27611

Purity Specification: 98%

Molecular Formula : C₄₇H₅₈N₁₂O₆

CASNumber: 140703-51-1

FormulaWeight: 887.040


Storage: 2-8

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816