

Specification Sheet



Product Name: 5,6-trans-Vitamin D3

Item No.: MED27794

Purity Specification: 98%

Molecular Formula : C₂₇H₄₄O

CASNumber: 22350-41-0

FormulaWeight: 384.638


Storage: -20

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816