

$^1\text{H}$  NMR (400 MHz,  $\text{DMSO}-d_6$ )  $\delta$  ppm 0.87 - 1.66 (m, 30 H) 2.85 (ddd,  $J=13.81, 6.87, 6.81$  Hz, 1 H) 2.91 - 3.08 (m, 2 H) 3.38 - 3.43 (m, 2 H) 3.83 (s, 2 H) 7.00 (s, 2 H) 7.16 - 7.49 (m, 3 H) 12.95 (br. s., 1 H)

