

Specification Sheet



Product Name: TW-37

Item No.: MED10109

Purity Specification: >96%


Molecular Formula : C₃₃H₃₅NO₆S

CASNumber: 877877-35-5 **FormulaWeight:** 573.699

Storage: -20

Test	Specification	
	Min.	Max.
PURITY	96%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816