

Specification Sheet



Product Name: Danoprevir

Item No.: MED10288

Purity Specification: >99%

Molecular Formula : C₃₅H₄₆FN₅O₉S

CASNumber: 850876-88-9

FormulaWeight: 731.831


Storage: -20°C

Test

Specification

	Min.	Max.
PURITY	99%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816