

Specification Sheet



Product Name: ML-323

Item No.: MED10598

Purity Specification: >99%


Molecular Formula : C₂₃H₂₄N₆

CASNumber: 1572414-83-5 **FormulaWeight:** 384.477

Storage: -20

Test	Specification	
	Min.	Max.
PURITY	99%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816