

Specification Sheet



Product Name: Tapentadol

Item No.: MED10893

Purity Specification: >98%

Molecular Formula : C₁₄H₂₄CLNO

CASNumber: 175591-09-0 **FormulaWeight:** 257.799


Storage: -20°C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816