

# Specification Sheet



**Product Name:** HMP Linker

*Item No.:* MED11017

**Purity Specification:** 98%

**Molecular Formula :** C<sub>9</sub>H<sub>10</sub>O<sub>4</sub>

**CASNumber:** 68858-21-9

**FormulaWeight:** 182.173


**Storage:** -20 ° C

## Test

## Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: [www.med-bio.cn](http://www.med-bio.cn) Tel: 400-186-1816