

Specification Sheet



Product Name: DAPK Substrate Peptide

Item No.: MED11243

Purity Specification: $\geq 98\%$

Molecular Formula :

CASNumber: 386769-53-5 *FormulaWeight:*

Storage: Desiccate at -20 °C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816