

Specification Sheet



Product Name: PhiKan 083

Item No.: MED11252

Purity Specification: $\geq 98\%$


Molecular Formula : C₁₆H₁₉CLN₂

CASNumber: 880813-36-5 **FormulaWeight:** 274.788

Storage: Store at 4 ° C

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816