

# Specification Sheet



**Product Name:** Lenalidomide hemihydrate

*Item No.:* MED11277

**Purity Specification:** 98%

**Molecular Formula :** C<sub>26</sub>H<sub>28</sub>N<sub>6</sub>O<sub>7</sub>


**CASNumber:** 847871-99-2

**FormulaWeight:** 536.537

**Storage:** -20 ° C

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: [www.med-bio.cn](http://www.med-bio.cn) Tel: 400-186-1816