

Specification Sheet



Product Name: PD-1/PD-L1 Inhibitor 3

Item No.: MED11304

Purity Specification: 95%

Molecular Formula :

CASNumber: 1629654-95-0 **FormulaWeight:**

Storage: -20 ° C

Test	Specification	
	Min.	Max.
PURITY	95%	100%

Structure :

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816