

Specification Sheet



Product Name: Angiotensin III (human, mouse)

Item No.: MED11430

Purity Specification: 98%

Molecular Formula : C₄₆H₆₆N₁₂O₉

CAS Number: 13602-53-4

Formula Weight: 931.091


Storage: Desiccate at -20 °C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816