

Specification Sheet



Product Name: Senexin B

Item No.: MED11494

Purity Specification: $\geq 98\%$

Molecular Formula : H₃N

CASNumber: 1449228-40-3 *FormulaWeight:* 17.031

Storage: -20 ° C

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816