

Specification Sheet



Product Name: ENMD-2076 L-(+)-Tartaric acid

Item No.: MED11671

Purity Specification: $\geq 98\%$

Molecular Formula : C₂₅H₃₁N₇O₆

CASNumber: 1291074-87-7 **FormulaWeight:** 525.557


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816