

Specification Sheet



Product Name: Calcineurin Autoinhibitory Peptide

Item No.: MED11719

Purity Specification: 98%


Molecular Formula : C₁₂₄H₂₀₅N₃₉O₃₉S₂

CASNumber: 148067-21-4 **FormulaWeight:** 2930.322

Storage: Desiccate at -20 °C

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816